



SCITECH®

Innovation for Life



V A S C U L A R

SNEAK
MICROCATHETER



Characteristics:

- Hydrophilic coating to ensure smooth introduction to navigate through tortuous anatomy
- Coil braided design provides outstanding pushability and kink resistance
- Radiopaque marker bands

Sneak Microcatheter Compatibility

Ref.	Catheter	Size	Marker Bands	Total Length	ID	Embsoft Range (µm)				
						100-300	300-500	500-700	700-900	900-1100
112065	SNEAK 28	2,6F	1	130cm	0,71mm/0,028"	●	●	●	●	
112063	SNEAK 24	2,4F	2	150cm	0,61mm/0,024"	●	●	●		
112059	SNEAK 21	2,1F	2	150cm	0,53mm/0,021"	●	●	●		
112054	SNEAK 17	1,9F	2	150cm	0,43mm/0,017"	●	●			